ORBIS™ MICROBUBBLE GENERATOR

INSTRUCTIONS FOR USE

DEVICE DESCRIPTION

The Orbis Microbubble Generator is designed solely for parenteral administration of air microbubble contrast via Luer Lock attachment to an intravenous line or needle and to a 10mL prefilled saline syringe.

It is recommended to inject a maximum of 4 times within 1 hour with the Orbis Microbubble Generator with 5 minutes recommended between injections. The device is provided sterile.

TABLE 1. Compatible Connections

20-22G needles
20-22G intravenous lines

INDICATIONS FOR USE

The Orbis Microbubble Generator is intended to be used for the intra-venous administration of contrast media via an intravenous catheter or needle. The Orbis Microbubble Generator is intended to attach to a pre-loaded saline syringe with contrast media and is used during ultrasonography procedures.

ORBIS MICROBUBBLE GENERATOR


FIG. 1: ORBIS MICROBUBBLE GENERATOR

DEVICE PARTS (FIG. 1)

1. Luer Lock Prefilled Syringe Connector
2. Intravenous Catheter or Needle Connector

HOW SUPPLIED

Sterile for single use only

STORAGE

Store at controlled room temperature: 59-86°F (15-30°C)

CONTRAINDICATIONS

- No known contraindications

WARNINGS

- Single patient use only. Reuse may cause an infection.
- Do not use if sterile package is opened or damaged. Do not re-sterilize.
- Do not use if device is damaged.
- Do not use if the device has been dropped, as the device may be broken even if a break cannot be seen.
- Do not use the device after the expiration date on the label.
- Visually inspect and do not use if there is any visible particulate matter.
- Ensure use of aseptic techniques throughout the procedure.
- Use only with compatible (20-22G) needles or directly into intravenous lines.

INSTRUCTIONS FOR USE

1. Inspect the package prior to use.
2. Verify that you have correct intravenous (IV) access and prepare the IV hub for connection to the device OR prepare a compatible 20-22G needle.
3. Peel open package and remove the device.
4. Prepare a 10mL prefilled saline syringe per the manufacturer's recommendations. Connect the Orbis to the 10mL prefilled saline syringe.
5. Connect the Orbis to the port of an IV line or needle of 20-22G, taking care that there is no touch contamination of the connection.
6. Using constant pressure at an angle of approximately 30 degrees, push the plunger rod as fast as possible and deliver contrast until the plunger rod bottoms out.
7. After use, disconnect the Orbis and dispose of product and packaging in accordance with hospital, administrative and/or local government policy. **DO NOT REUSE.**
8. Confirmation of full clearance of contrast is recommended prior to repeating an injection.

PACKAGING, SHELF-LIFE AND STORAGE

The device packaging is designed to minimize exposure to the atmosphere and to provide for aseptic product transfer. Store the product in the unopened package until time of use. Contents are sterile prior to the expiration date of the package if the package is undamaged and unopened. Do not re-sterilize. The expiration date is marked on the carton box and on the individual devices. The package should be stored in a cool, dry location. Store at 59-86°F (15-30°C). Dispose of used product and packaging in accordance with local, state and federal regulations.

SYMBOLS GLOSSARY

	Model Number		Do not use if package is opened or damaged		Use By Date
	Lot Number		Do Not Reuse		Sterilized Using Irradiation
	Non-Pyrogenic		Prescription Only		Temperature Range between 59-86°F (15-30°C)
	Consult Instructions for Use		Warning		Not made with natural rubber latex
	Do Not Re-sterilize		Manufactured By		Quantity

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