Quick Reference Guide

INDICATIONS FOR USE

The Orbis Microbubble Generator is intended to be used for the intravenous administration of contrast media via an intravenous catheter or needle. The Orbis Microbubble Generator is intended to attach to a pre-loaded saline syringe with contrast media and is used during ultrasonography procedures. **Please refer to the IFU for full instructions, contraindications and warnings.**

INSTRUCTIONS FOR USE

- 1 Verify that you have correct intravenous (IV) access and prepare the IV hub for connection to the device OR prepare a compatible 20-22G needle.
- 2 After inspecting package, peel open and remove the device using sterile technique.
- 3 Prepare a 10mL prefilled saline syringe per the manufacturer's recommendations. Connect the Orbis to the 10mL prefilled saline syringe [Fig 1].
- 4 Connect the Orbis to the port of an IV line or needle of 20-22G, taking care that there is no touch contamination of the connection.
- 5 Using constant pressure at an angle of approximately 30 degrees, push the plunger rod as fast as possible and deliver contrast until the plunger rod bottoms out [Fig 2.].
- 6 After use, disconnect the Orbis and dispose of product and packaging in accordance with hospital, administrative and/or local government policy. **DO NOT REUSE.**
- 7 Confirmation of full clearance of contrast is recommended prior to repeating an injection.

Fig 1. Connection of Orbis to saline syringe.


Fig 2. Hold syringe at 30 degree angle and push plunger to deliver contrast.

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